

TuneMyBG LLM reproducibility — plain-English summary

Package generated 2026-08-21 covering 2026-08-07 to 2026-08-21 (14 days). Glucose: mean 121.6 mg/dL, time in range 86.0%, below range 5.1%, CV 32.6%. Carb entries logged: 0.

Each run is a brand-new conversation: package attached with Prompt 1, then Prompts 2, 3 and 4 in turn, then the app’s correction prompt once if the JSON failed the app’s checks. AAPS defaults are from nightscout/AndroidAPS master @ 598e2eb (2026-08-02).

gemini-3.6-flash — 50 runs (real_gemini36flash_50)

At a glance

- App outcome over 50 completed runs: accepted first time 50; accepted after the repair prompt 0; invalid JSON with no prompt offered 0; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 0.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 46/50 (92%) of runs first time.
- Settings changed per run: 1–6 of 23 (most runs: 3).
- 26 settings were left unchanged in every single run; 10 were changed in at least one run, of which 4 were changed in more than half the runs.
- Headline profile recommendation focus: target in 39/50 (78%), basal in 11/50 (22%).
- Profile recommendation decision: change ×44, keep ×6.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	49/50 (98%)	4–15 U; most often 7 U (27 of 49); 7 distinct values	lower
AAPS / max basal and max IOB limits / max_basal_u_per_hour	12 U/h	1.0 U/h (Simple Mode recalculates this from the profile)	48/50 (96%)	2.5–4 U/h; most often 3.5 U/h (19 of 48); 5 distinct values	lower
Profile / target 00:00	90 mg/dL	—	39/50 (78%)	always 100 mg/dL	higher
Profile / target 23:00	90 mg/dL	—	39/50 (78%)	always 100 mg/dL	higher
AAPS / Autosens and Dynamic ISF / dynamic_isf_adjustment_factor_percent	70 %	100 %	13/50 (26%)	always 50 %	lower
Profile / basal 00:00	0.6 U/h	—	11/50 (22%)	always 0.5 U/h	lower
AAPS / SMB strength and frequency / max_uam_smb_basal_minutes	20 min	30 min	4/50 (8%)	always 15 min	lower
Profile / basal 19:00	0.7 U/h	—	3/50 (6%)	0.5–0.6 U/h; most often 0.5 U/h (2 of 3); 2 distinct values	lower
AAPS / carbohydrate absorption model / min_5m_carbimpact	8 mg/dL/5 min	8 mg/dL/5 min	1/50 (2%)	always 4 mg/dL/5 min	lower
Profile / basal 18:00	0.6 U/h	—	1/50 (2%)	always 0.5 U/h	lower

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 08:00 100 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false
- **AAPS carb_absorption:** meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: smb x22, target x16, safety_limits x5, basal x3, aaps x3, safety x1.

- Reduce AAPS Max IOB and Max Basal Safety Limits x9
- Raise Overnight Target to 100 mg/dL x7
- Reduce Max IOB and Max Basal Safety Limits x5
- Tighten AAPS Max IOB and Max Basal Safety Limits x4
- Raise Overnight Target to Reduce Nocturnal Hypoglycemia x2
- Reduce Overnight Basal Rate x2
- Raise Overnight Glucose Target x2
- Reduce AAPS Automation Safety Caps x1
- ... and 18 other titles

gemini-3.1-pro-preview — 50 runs (real_gemini31pro_50)

At a glance

- App outcome over 50 completed runs: accepted first time 50; accepted after the repair prompt 0; invalid JSON with no prompt offered 0; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 0.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 47/50 (94%) of runs first time.
- Settings changed per run: 1–4 of 23 (most runs: 3).
- 29 settings were left unchanged in every single run; 7 were changed in at least one run, of which 4 were changed in more than half the runs.
- Headline profile recommendation focus: target in 46/50 (92%), basal in 4/50 (8%).
- Profile recommendation decision: change x45, keep x5.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	45/50 (90%)	3–10 U; most often 7 U (16 of 45); 7 distinct values	lower
Profile / target 00:00	90 mg/dL	—	45/50 (90%)	100–110 mg/dL; most often 100 mg/dL (24 of 45); 2 distinct values	higher

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / target 23:00	90 mg/dL	—	45/50 (90%)	100–110 mg/dL; most often 100 mg/dL (24 of 45); 2 distinct values	higher
AAPS / max basal and max IOB limits / max_basal_u_per_hour	12 U/h	1.0 U/h (Simple Mode recalculates this from the profile)	41/50 (82%)	2.5–4 U/h; most often 3 U/h (29 of 41); 5 distinct values	lower
AAPS / SMB strength and frequency / smb_interval_minutes	3 min	3 min	3/50 (6%)	5–15 min; most often 15 min (2 of 3); 2 distinct values	higher
Profile / basal 00:00	0.6 U/h	—	1/50 (2%)	always 0.5 U/h	lower
Profile / target 08:00	100 mg/dL	—	1/50 (2%)	always 110 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: target ×17, automation ×9, safety_limits ×8, aaps.core.safety_limits ×6, aaps ×4, safety ×3, automation_safety ×1, automation_sensitivity_limits ×1, limits ×1.

- Reduce AAPS Safety Limits ×3
- Reduce Max IOB and Max Basal Limits ×3
- Reduce Max IOB and Max Basal limits ×2
- Raise Overnight Target for Safety ×2
- Raise Overnight Glucose Targets ×2
- Restrict Max IOB and Max Basal ×2
- Create a safety buffer against overnight lows ×1
- Raise overnight glucose target for safety ×1
- ... and 34 other titles

gemini-3.5-flash-lite — 49 runs (real_gemini35flashLite_50)

At a glance

- App outcome over 50 completed runs: accepted first time 38; accepted after the repair prompt 1; invalid JSON with no prompt offered 10; schema stop with no prompt offered 1; still rejected after repair 0. Dead ends (no way forward in the app): 11.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 32/49 (65%) of runs first time.
- Settings changed per run: 0–2 of 23 (most runs: 0).
- 33 settings were left unchanged in every single run; 3 were changed in at least one run, of which 0 were changed in more than half the runs.

- Headline profile recommendation focus: basal in 41/49 (84%), isf in 3/49 (6%), dia in 2/49 (4%), cr in 2/49 (4%), target in 1/49 (2%).
- Profile recommendation decision: keep ×43, change ×6.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	8/49 (16%)	6–15 U; most often 6 U (5 of 8); 4 distinct values	lower
Profile / basal 00:00	0.6 U/h	—	5/49 (10%)	always 0.5 U/h	lower
Profile / target 00:00	90 mg/dL	—	1/49 (2%)	always 95 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 08:00 100 mg/dL; 23:00 90 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: basal ×22, meal ×7, cr ×6, smb ×5, aaps ×5, meal_strategy ×1, profile ×1, target ×1, safety ×1.

- Review Max IOB Safety Limit ×4
- Verify Overnight Basal Rates ×3
- Initiate Carbohydrate Logging ×3
- Introduce Manual Carbohydrate Logging ×3
- Verify Overnight Basal and Fasting Stability ×3
- Reduce Overnight Basal Rates ×2
- Reduce Max IOB Safety Limit ×2
- Reduce Max IOB Limit for Safety ×2
- ... and 27 other titles

gpt-5.6-sol — 50 runs (real_gpt56sol_med_50)

Reasoning effort: medium.

At a glance

- App outcome over 50 completed runs: accepted first time 50; accepted after the repair prompt 0; invalid JSON with no prompt offered 0; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 0.
- Strict compliance with Prompt 4's own rules (which the app does not fully check): 50/50 (100%) of runs first time.
- Settings changed per run: 0–2 of 23 (most runs: 0).
- 34 settings were left unchanged in every single run; 2 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: target in 39/50 (78%), basal in 11/50 (22%).
- Profile recommendation decision: keep ×32, change ×18.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / target 00:00	90 mg/dL	—	18/50 (36%)	always 100 mg/dL	higher
Profile / target 23:00	90 mg/dL	—	2/50 (4%)	always 100 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 00:00 0.6 U/h; 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 08:00 100 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h; max_job_units = 25 U
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: target ×28, basal ×9, profile ×8, safety_and_data ×2, data_and_profile ×2, safety ×1.

- Primary: review a higher 00:00 target ×3
- Keep the profile while verifying the overnight low pattern ×2
- Verify the overnight low pattern before changing the profile ×2
- Primary: keep the profile stable and verify the low pattern first ×1
- Review the 00:00 glucose target first ×1
- Primary: verify the context of recurrent lows before changing the profile ×1
- Primary: review the midnight target ×1
- Keep the profile and verify the overnight pattern first ×1
- ... and 38 other titles

Reasoning effort: medium.

At a glance

- App outcome over 50 completed runs: accepted first time 48; accepted after the repair prompt 1; invalid JSON with no prompt offered 1; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 1.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 23/50 (46%) of runs first time.
- Settings changed per run: 0–1 of 23 (most runs: 0).
- 35 settings were left unchanged in every single run; 1 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 49/50 (98%), cr in 1/50 (2%).
- Profile recommendation decision: keep ×48, change ×2.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / basal 00:00	0.6 U/h	—	2/50 (4%)	always 0.55 U/h	lower

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 00:00 90 mg/dL; 08:00 100 mg/dL; 23:00 90 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h; max_job_units = 25 U
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: basal ×34, profile ×13, basal_and_targets ×1, meal ×1, profile_and_data ×1.

- Keep the current profile unchanged and verify with fasting data first ×1
- Keep the current profile unchanged and verify with separated no-meal observations ×1
- Keep profile unchanged pending one fasting verification ×1
- Keep the current profile and verify with one focused test ×1
- Keep the current profile and verify with meal-linked observation first ×1
- Keep the current profile unchanged and verify with one clean fasting-like block ×1
- Keep the current profile unchanged and verify with a meal-logged day before any basal, CR, ISF, target, or DIA revision ×1
- Keep the current profile unchanged for now ×1
- ... and 42 other titles

claude-opus-5 — 19 runs (real_opus5)

At a glance

- App outcome over 20 completed runs: accepted first time 10; accepted after the repair prompt 1; invalid JSON with no prompt offered 9; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 9.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 5/19 (26%) of runs first time.
- Settings changed per run: 1–3 of 23 (most runs: 1).
- 31 settings were left unchanged in every single run; 5 were changed in at least one run, of which 1 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 13/19 (68%), dia in 3/19 (16%), isf in 3/19 (16%).
- Profile recommendation decision: keep x19.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	15/19 (79%)	6–10 U; most often 8 U (8 of 15); 4 distinct values; a reachable backstop with clear headroom above the observed peak of 5.94 u, for example 8-10 u, to be confirmed with the care team x1	lower
AAPS / Autosens and Dynamic ISF / dynamic_isf_adjustment_factor_perce nt	70 %	100 %	6/19 (32%)	55–60 %; most often 60 % (5 of 6); 2 distinct values	lower
AAPS / SMB strength and frequency / max_uam_smb_basal_minutes	20 min	30 min	4/19 (21%)	always 15 min	lower
AAPS / max basal and max IOB limits / max_basal_u_per_hour	12 U/h	1.0 U/h (Simple Mode recalculates this from the profile)	2/19 (11%)	4–5 U/h; most often 4 U/h (1 of 2); 2 distinct values	lower
AAPS / SMB strength and frequency / smb_interval_minutes	3 min	3 min	1/19 (5%)	always 5 min	higher

Kept unchanged in every run

- **Basal rates:** 00:00 0.6 U/h; 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 00:00 90 mg/dL; 08:00 100 mg/dL; 23:00 90 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** max_smb_basal_minutes = 15 min
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: smb x14, isf x3, safety_limits x1, meal x1.

- Primary: review the max IOB ceiling as the first and only change this cycle x1
- Reduce the maximum insulin-on-board ceiling as the single first change x1
- Primary: discuss reducing the unannounced-meal microbolus ceiling as the single first change x1
- Discuss reducing max IOB as the first, independent step x1
- Reduce the maximum insulin-on-board ceiling to a value that can actually act x1
- Discuss tightening the maximum insulin on board ceiling as the single primary change x1
- Primary: reduce the dynamic ISF adjustment factor as the single first change x1
- Primary: discuss aligning the UAM microbolus basal-minutes cap with the non-UAM cap x1
- ... and 11 other titles

claude-sonnet-5 — 20 runs (real_sonnet5_20)

At a glance

- App outcome over 20 completed runs: accepted first time 18; accepted after the repair prompt 1; invalid JSON with no prompt offered 1; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 1.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 16/20 (80%) of runs first time.
- Settings changed per run: 0–1 of 23 (most runs: 0).
- 35 settings were left unchanged in every single run; 1 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 14/20 (70%), target in 3/20 (15%), dia in 2/20 (10%), cr in 1/20 (5%).
- Profile recommendation decision: keep x19, change x1.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / basal 16:00	0.5 U/h	—	1/20 (5%)	approximately 0.6 u/h (modest increase, verify stepwise) x1	—

Kept unchanged in every run

- **Basal rates:** 00:00 0.6 U/h; 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 00:00 90 mg/dL; 08:00 100 mg/dL; 23:00 90 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h; max_iob_units = 25 U
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %

- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: meal ×6, safety ×3, data ×2, basal ×2, dia ×2, target ×1, cr ×1, data_quality ×1, safety_job ×1, data_logging ×1.

- Prioritize understanding recurrent overnight hypoglycemia before any other change ×1
- Begin logging carbohydrate and e-carb entries ×1
- Begin logging carbohydrate intake to resolve competing explanations ×1
- Establish basic carbohydrate logging before further profile or automation changes ×1
- Clarify the effective maximum basal ceiling before any basal change ×1
- Establish carbohydrate logging before further profile or automation changes ×1
- Investigate the recurring correction-then-delayed-hypoglycemia pattern before any other change ×1
- Establish basic carbohydrate logging at recurring meal windows ×1
- ... and 12 other titles

claude-haiku-4-5 — 27 runs (real_haiku45_50)

At a glance

- App outcome over 28 completed runs: accepted first time 13; accepted after the repair prompt 0; invalid JSON with no prompt offered 2; schema stop with no prompt offered 0; still rejected after repair 13. Dead ends (no way forward in the app): 2.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 1/27 (4%) of runs first time.
- Settings changed per run: 0–5 of 23 (most runs: 1).
- 23 settings were left unchanged in every single run; 13 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 12/27 (44%), ISF in 1/27 (4%), basal_cr_isf_dia_validation_sequence in 1/27 (4%), basal, target, CR, ISF, DIA in 1/27 (4%), carb logging, then isf verification in 1/27 (4%), carbohydrate_ratio_and_meal_strategy in 1/27 (4%), insulin_duration in 1/27 (4%), carbohydrate_logging_and_meal_bolus_execution in 1/27 (4%), DIA and basal timing (safety priority), conditional midday basal increase (medium priority, contingent on ISF verification) in 1/27 (4%), basal_and_targets_and_cr_and_isf_and_dia in 1/27 (4%), carb_ratio in 1/27 (4%), isf in 1/27 (4%), meal_strategy_and_carb_logging in 1/27 (4%), basal, cr, isf, dia, target in 1/27 (4%), meal_bolus_timing_and_carbohydrate_entry in 1/27 (4%), Profile defects: flat ISF 56 mg/dL/U (too high, masked by Dynamic ISF compression to 2–4); overnight basal 0.6 U/h (excessive, causing hour 1 low 90.6 mg/dL, 18.5% below range); afternoon CR 13.2 g/U (loosest segment, coinciding with sustained hyperglycemia 130–162 mg/dL, 26% above range). GLP-1 post-2026-08-14 extends meal absorption to 240–360m, but CR and ISF unchanged; reactive bolusing (no prebolus, no carb logging) allows meal boluses to peak at 60m while carbs continue absorbing, causing rebounds to 188–235 mg/dL. AAPS automation (SMB, UAM, temp basal, Dynamic ISF, safety limits) is appropriately configured; no automation defects. Automation is compensating for profile-level errors in real-time. in 1/27 (4%).
- Profile recommendation decision: change ×17, keep ×10.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / basal 10:00	0.7 U/h	—	10/27 (37%)	0.8–0.85 U/h; most often 0.8 U/h (6 of 10); 2 distinct values	higher
Profile / basal 16:00	0.5 U/h	—	6/27 (22%)	0.6–0.7 U/h; most often 0.7 U/h (3 of 6); 3 distinct values	higher
Profile / basal 00:00	0.6 U/h	—	3/27 (11%)	0.5–0.55 U/h; most often 0.5 U/h (2 of 3); 2 distinct values	lower

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / DIA	10 h	— (AAPS defaultDIA = 5.0 h (Constants.kt))	3/27 (11%)	5.5–6 h; most often 5.5 h (2 of 3); 2 distinct values	lower
Profile / ISF 00:00	56 mg/dL/U	—	3/27 (11%)	40–110 mg/dL/U; most often 90 mg/dL/U (1 of 3); 3 distinct values	mixed
AAPS / carbohydrate absorption model / min_5m_carbimpact	8 mg/dL/5 min	8 mg/dL/5 min	2/27 (7%)	5–12 mg/dL/5 min; most often 12 mg/dL/5 min (1 of 2); 2 distinct values	mixed
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	2/27 (7%)	8–20 U; most often 8 U (1 of 2); 2 distinct values	lower
AAPS / Autosens and Dynamic ISF / dynamic_isf_adjustment_factor_percent	70 %	100 %	1/27 (4%)	always 85 %	higher
Profile / basal 19:00	0.7 U/h	—	1/27 (4%)	always 0.6 U/h	lower
Profile / CR 08:00	13.2 g/U	—	1/27 (4%)	always 11.5 g/U	lower
Profile / CR 16:00	10.5 g/U	—	1/27 (4%)	always 9.5 g/U	lower
Profile / target 00:00	90 mg/dL	—	1/27 (4%)	always 100 mg/dL	higher
Profile / target 23:00	90 mg/dL	—	1/27 (4%)	always 100 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 06:00 0.8 U/h; 18:00 0.6 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 20:00 12.6 g/U
- **Targets:** 08:00 100 mg/dL
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false
- **AAPS carb_absorption:** meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: basal x9, meal_strategy x4, ISF x1, Data and meal strategy x1, safety x1, meal strategy x1, Carbohydrate logging (highest-leverage behavior change) x1, meal_logging x1, dia x1, behavioral_and_logging x1, Data Integrity – ISF Verification x1, carb_logging x1, isf x1, meal_bolus_strategy x1, meal_carbs x1, profile x1.

- Correct ISF underestimation from 56 to 90 mg/dL/U x1
- Restart carb logging immediately for 3 complete days x1
- Reduce max_iob from 25 U to 6-8 U immediately x1
- Enable carb logging (Nightscout treatment entries) immediately for all meals and snacks; log carb amount (grams) and meal time. Continue for minimum 7 days to establish baseline. x1
- Reduce basal 00:00 from 0.6 U/h to 0.55 U/h or 0.5 U/h to address hour 1 lows (18.5% TBR, nadir 42–55 mg/dL). x1
- Raise basal in 10:00–16:00 afternoon window to address persistent hyperglycemia. x1
- Implement meal carbohydrate logging x1
- Reduce overnight and evening basal segments to address recurrent hypoglycemia x1
- ... and 19 other titles

deepseek/deepseek-v4-pro — 50 runs (real_deepseek_v4pro_50)

At a glance

- App outcome over 50 completed runs: accepted first time 45; accepted after the repair prompt 3; invalid JSON with no prompt offered 0; schema stop with no prompt offered 2; still rejected after repair 0. Dead ends (no way forward in the app): 2.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 33/50 (66%) of runs first time.
- Settings changed per run: 0–6 of 23 (most runs: 2).
- 18 settings were left unchanged in every single run; 18 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 28/50 (56%), isf in 17/50 (34%), dia in 4/50 (8%), target in 1/50 (2%).
- Profile recommendation decision: change ×39, keep ×11.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	22/50 (44%)	2–12 U; most often 5 U (7 of 22); 7 distinct values	lower
Profile / basal 00:00	0.6 U/h	—	22/50 (44%)	0.5–0.55 U/h; most often 0.5 U/h (20 of 22); 2 distinct values	lower
Profile / ISF 00:00	56 mg/dL/U	—	20/50 (40%)	30–180 mg/dL/U; most often 70 mg/dL/U (5 of 20); 14 distinct values	mixed
AAPS / max basal and max IOB limits / max_basal_u_per_hour	12 U/h	1.0 U/h (Simple Mode recalculates this from the profile)	9/50 (18%)	2–4 U/h; most often 4 U/h (4 of 9); 3 distinct values	lower
AAPS / Autosens and Dynamic ISF / dynamic_isf_adjustment_factor_percent	70 %	100 %	9/50 (18%)	0–100 %; most often 100 % (4 of 9); 5 distinct values	mixed
Profile / DIA	10 h	— (AAPS defaultDIA = 5.0 h (Constants.kt))	8/50 (16%)	5–8 h; most often 6 h (4 of 8); 4 distinct values	lower
AAPS / SMB strength and frequency / max_uam_smb_basal_minutes	20 min	30 min	5/50 (10%)	10–15 min; most often 15 min (3 of 5); 2 distinct values	lower
AAPS / carbohydrate absorption model / min_5m_carbimpact	8 mg/dL/5 min	8 mg/dL/5 min	3/50 (6%)	3–5 mg/dL/5 min; most often 5 mg/dL/5 min (2 of 3); 2 distinct values	lower
Profile / basal 19:00	0.7 U/h	—	2/50 (4%)	always 0.6 U/h	lower
AAPS / carbohydrate absorption model / meal_max_absorption_hours	7 h	6 h	1/50 (2%)	always 5 h	lower
AAPS / SMB activation conditions / smb_always_enabled	true	true	1/50 (2%)	false ×1	—
AAPS / SMB activation conditions / smb_with_temp_target_enabled	true	true	1/50 (2%)	false ×1	—
AAPS / SMB strength and frequency / smb_interval_minutes	3 min	3 min	1/50 (2%)	always 10 min	higher

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / SMB strength and frequency / max_smb_basal_minutes	15 min	30 min	1/50 (2%)	always 10 min	lower
AAPS / SMB and UAM / uam_enabled	true	true	1/50 (2%)	false x1	—
Profile / basal 05:00	0.7 U/h	—	1/50 (2%)	always 0.6 U/h	lower
Profile / target 00:00	90 mg/dL	—	1/50 (2%)	always 100 mg/dL	higher
Profile / target 23:00	90 mg/dL	—	1/50 (2%)	always 100 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Targets:** 08:00 100 mg/dL
- **AAPS smb_uam:** smb_enabled = true
- **AAPS smb_activation:** smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_high_temp_target_enabled = false
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false

Primary recommendation (first item the model listed)

By area: basal x18, isf x16, aaps_safety_limits x2, dia x2, meal_strategy x1, data_foundation x1, aaps_safety x1, aaps x1, automation_sensitivity_limits x1, uam_aggressiveness x1, UAM SMB delivery x1, automation_sensitivity x1, meal x1, target x1, aaps_uam x1.

- Reduce overnight basal rate x6
- Increase Insulin Sensitivity Factor (ISF) x2
- Reduce overnight basal rate 00:00-05:00 x2
- Begin carbohydrate entries for all meals x1
- Reduce overnight basal (00:00–05:00) x1
- Record all carbohydrate intake consistently x1
- Correct Insulin Sensitivity Factor (ISF) x1
- Reduce overnight basal rate (00:00–05:00) x1
- ... and 34 other titles

x-ai/grok-4.6 — 50 runs (real_grok46_50)

At a glance

- App outcome over 50 completed runs: accepted first time 50; accepted after the repair prompt 0; invalid JSON with no prompt offered 0; schema stop with no prompt offered 0; still rejected after repair 0. Dead ends (no way forward in the app): 0.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 49/50 (98%) of runs first time.
- Settings changed per run: 0–2 of 23 (most runs: 0).
- 34 settings were left unchanged in every single run; 2 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: target in 33/50 (66%), basal in 17/50 (34%).
- Profile recommendation decision: keep x40, change x10.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / target 00:00	90 mg/dL	—	10/50 (20%)	always 100 mg/dL	higher
Profile / target 23:00	90 mg/dL	—	6/50 (12%)	always 100 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 00:00 0.6 U/h; 05:00 0.7 U/h; 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Insulin sensitivity:** 00:00 56 mg/dL/U
- **Targets:** 08:00 100 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h; max_job_units = 25 U
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: target x33, basal x16, profile x1.

- Keep the profile and first separate overnight lows from leftover IOB x2
- Keep the profile and do not add insulin for afternoon highs x1
- Review raising only the 00:00 target from 90 to 100 mg/dL x1
- Keep the profile and first separate overnight IOB from target and basal x1
- Keep the profile and separate leftover IOB from overnight landing x1
- Keep the current profile and verify nights without leftover meal IOB x1
- Keep the overnight profile and compare midnight IOB with overnight lows x1
- Keep the profile and verify overnight lows against leftover IOB x1
- ... and 41 other titles

meta-llama/llama-4-maverick — 49 runs (real_llama4maverick_50)

At a glance

- App outcome over 50 completed runs: accepted first time 19; accepted after the repair prompt 25; invalid JSON with no prompt offered 2; schema stop with no prompt offered 0; still rejected after repair 4. Dead ends (no way forward in the app): 2.
- Strict compliance with Prompt 4’s own rules (which the app does not fully check): 19/49 (39%) of runs first time.
- Settings changed per run: 0–4 of 23 (most runs: 1).
- 30 settings were left unchanged in every single run; 6 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 43/49 (88%), isf in 4/49 (8%), cr in 1/49 (2%), aaps in 1/49 (2%).

- Profile recommendation decision: change x29, keep x20.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / basal 00:00	0.6 U/h	—	19/49 (39%)	0.55–0.65 U/h; most often 0.65 U/h (15 of 19); 2 distinct values	mixed
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	18/49 (37%)	15–20 U; most often 15 U (12 of 18); 2 distinct values	lower
Profile / ISF 00:00	56 mg/dL/U	—	4/49 (8%)	50–55 mg/dL/U; most often 50 mg/dL/U (3 of 4); 2 distinct values	lower
Profile / basal 16:00	0.5 U/h	—	2/49 (4%)	0.55–0.6 U/h; most often 0.55 U/h (1 of 2); 2 distinct values	higher
Profile / basal 06:00	0.8 U/h	—	1/49 (2%)	always 0.75 U/h	lower
Profile / CR 00:00	13 g/U	—	1/49 (2%)	always 12.5 g/U	lower

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 10:00 0.7 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Targets:** 00:00 90 mg/dL; 08:00 100 mg/dL; 23:00 90 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min
- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: basal x26, AAPS safety limits x7, profile x5, AAPS x4, meal strategy x3, Meal Logging x1, meal logging x1, Meal Strategy x1.

- Adjust basal rates x6
- Verify basal rates x4
- Adjust max IOB x4
- Adjust max IOB limit x3
- Review basal rates x3
- Verify and adjust basal rates x2
- Start logging meal data x2
- Reduce max IOB x2
- ... and 20 other titles

Model comparison — how often each setting was changed

Setting	Current	AAPS default	gemini-3.6-flash	gemini-3.1-pro-preview	gemini-3.5-flash-lite	gpt-5.6-sol	gpt-5.4-mini	claude-opus-5	claude-sonnet-5	claude-haiku-4-5	deepseek/deepseek-v4-pro	x-ai/grok-4.6	meta-llama/llama-4-maverick
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U	98% (4–15)	90% (3–10)	16% (6–15)	never	never	79% (6–10)	never	7% (8–20)	44% (2–12)	never	37% (15–20)
AAPS / max basal and max IOB limits / max_basal_u_per_hour	12 U/h	1.0 U/h	96% (2.5–4)	82% (2.5–4)	never	never	never	11% (4–5)	never	never	18% (2–4)	never	never
Profile / target 00:00	90 mg/dL	—	78% (100)	90% (100–110)	2% (95)	36% (100)	never	never	never	4% (100)	2% (100)	20% (100)	never
Profile / target 23:00	90 mg/dL	—	78% (100)	90% (100–110)	never	4% (100)	never	never	never	4% (100)	2% (100)	12% (100)	never
AAPS / Autosens and Dynamic ISF / dynamic_isf_adjustment_factor_percent	70 %	100 %	26% (50)	never	never	never	never	32% (55–60)	never	4% (85)	18% (0–100)	never	never
Profile / basal 00:00	0.6 U/h	—	22% (0.5)	2% (0.5)	10% (0.5)	never	4% (0.55)	never	never	11% (0.5–0.55)	44% (0.5–0.55)	never	39% (0.55–0.65)
AAPS / SMB strength and frequency / max_uam_smb_basal_minutes	20 min	30 min	8% (15)	never	never	never	never	21% (15)	never	never	10% (10–15)	never	never
Profile / basal 19:00	0.7 U/h	—	6% (0.5–0.6)	never	never	never	never	never	never	4% (0.6)	4% (0.6)	never	never
AAPS / carbohydrate absorption model / min_5m_carb_impact	8 mg/dL/5 min	8 mg/dL/5 min	2% (4)	never	never	never	never	never	never	7% (5–12)	6% (3–5)	never	never
Profile / basal 18:00	0.6 U/h	—	2% (0.5)	never	never	never	never	never	never	never	never	never	never
AAPS / SMB strength and frequency / smb_interval_minutes	3 min	3 min	never	6% (5–15)	never	never	never	5% (5)	never	never	2% (10)	never	never

Setting	Current	AAPS default	gemini-3.6-flash	gemini-3.1-pro-preview	gemini-3.5-flash-lite	gpt-5.6-sol	gpt-5.4-mini	claude-opus-5	claude-sonnet-5	claude-haiku-4-5	deepseek/deepseek-v4-pro	x-ai/grok-4.6	meta-llama/llama-4-maverick
Profile / target 08:00	100 mg/dL	—	never	2% (110)	never	never	never	never	never	never	never	never	never
Profile / basal 16:00	0.5 U/h	—	never	never	never	never	never	never	5% (approximately 0.6 u/h (modest increase, verify stepwise))	22% (0.6–0.7)	never	never	4% (0.55–0.6)
Profile / basal 10:00	0.7 U/h	—	never	never	never	never	never	never	never	37% (0.8–0.85)	never	never	never
Profile / DIA	10 h	—	never	never	never	never	never	never	never	11% (5.5–6)	16% (5–8)	never	never
Profile / ISF 00:00	56 mg/dL/U	—	never	never	never	never	never	never	never	11% (40–110)	40% (30–180)	never	8% (50–55)
Profile / CR 08:00	13.2 g/U	—	never	never	never	never	never	never	never	4% (11.5)	never	never	never
Profile / CR 16:00	10.5 g/U	—	never	never	never	never	never	never	never	4% (9.5)	never	never	never
AAPS / carbohydrate absorption model / meal_max_absorption_hours	7 h	6 h	never	never	never	never	never	never	never	never	2% (5)	never	never
AAPS / SMB activation conditions / smb_always_enabled	true	true	never	never	never	never	never	never	never	never	2% (false)	never	never
AAPS / SMB activation conditions / smb_with_temp_target_enabled	true	true	never	never	never	never	never	never	never	never	2% (false)	never	never
AAPS / SMB strength and frequency / max_smb_basal_minutes	15 min	30 min	never	never	never	never	never	never	never	never	2% (10)	never	never
AAPS / SMB and UAM / uam_enabled	true	true	never	never	never	never	never	never	never	never	2% (false)	never	never
Profile / basal 05:00	0.7 U/h	—	never	never	never	never	never	never	never	never	2% (0.6)	never	never

[illegible]