

TuneMyBG LLM reproducibility — plain-English summary

Package generated 2026-08-21 covering 2026-08-07 to 2026-08-21 (14 days). Glucose: mean 121.6 mg/dL, time in range 86.0%, below range 5.1%, CV 32.6%. Carb entries logged: 0.

Each run is a brand-new conversation: package attached with Prompt 1, then Prompts 2, 3 and 4 in turn, then the app's correction prompt once if the JSON failed the app's checks. AAPS defaults are from nightscout/AndroidAPS master @ 598e2eb (2026-08-02).

deepseek/deepseek-v4-pro — 50 runs (real_deepseek_v4pro_50)

At a glance

- App outcome over 50 completed runs: accepted first time 45; accepted after the repair prompt 3; invalid JSON with no prompt offered 0; schema stop with no prompt offered 2; still rejected after repair 0. Dead ends (no way forward in the app): 2.
- Strict compliance with Prompt 4's own rules (which the app does not fully check): 33/50 (66%) of runs first time.
- Settings changed per run: 0–6 of 23 (most runs: 2).
- 18 settings were left unchanged in every single run; 18 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 28/50 (56%), isf in 17/50 (34%), dia in 4/50 (8%), target in 1/50 (2%).
- Profile recommendation decision: change ×39, keep ×11.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	22/50 (44%)	2–12 U; most often 5 U (7 of 22); 7 distinct values	lower
Profile / basal 00:00	0.6 U/h	—	22/50 (44%)	0.5–0.55 U/h; most often 0.5 U/h (20 of 22); 2 distinct values	lower
Profile / ISF 00:00	56 mg/dL/U	—	20/50 (40%)	30–180 mg/dL/U; most often 70 mg/dL/U (5 of 20); 14 distinct values	mixed
AAPS / max basal and max IOB limits / max_basal_u_per_hour	12 U/h	1.0 U/h (Simple Mode recalculates this from the profile)	9/50 (18%)	2–4 U/h; most often 4 U/h (4 of 9); 3 distinct values	lower
AAPS / Autosens and Dynamic ISF / dynamic_isf_adjustment_factor_percent	70 %	100 %	9/50 (18%)	0–100 %; most often 100 % (4 of 9); 5 distinct values	mixed
Profile / DIA	10 h	— (AAPS defaultDIA = 5.0 h (Constants.kt))	8/50 (16%)	5–8 h; most often 6 h (4 of 8); 4 distinct values	lower
AAPS / SMB strength and frequency / max_uam_smb_basal_minutes	20 min	30 min	5/50 (10%)	10–15 min; most often 15 min (3 of 5); 2 distinct values	lower
AAPS / carbohydrate absorption model / min_5m_carbimpact	8 mg/dL/5 min	8 mg/dL/5 min	3/50 (6%)	3–5 mg/dL/5 min; most often 5 mg/dL/5 min (2 of 3); 2 distinct values	lower
Profile / basal 19:00	0.7 U/h	—	2/50 (4%)	always 0.6 U/h	lower
AAPS / carbohydrate absorption model / meal_max_absorption_hours	7 h	6 h	1/50 (2%)	always 5 h	lower

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
AAPS / SMB activation conditions / smb_always_enabled	true	true	1/50 (2%)	false x1	—
AAPS / SMB activation conditions / smb_with_temp_target_enabled	true	true	1/50 (2%)	false x1	—
AAPS / SMB strength and frequency / smb_interval_minutes	3 min	3 min	1/50 (2%)	always 10 min	higher
AAPS / SMB strength and frequency / max_smb_basal_minutes	15 min	30 min	1/50 (2%)	always 10 min	lower
AAPS / SMB and UAM / uam_enabled	true	true	1/50 (2%)	false x1	—
Profile / basal 05:00	0.7 U/h	—	1/50 (2%)	always 0.6 U/h	lower
Profile / target 00:00	90 mg/dL	—	1/50 (2%)	always 100 mg/dL	higher
Profile / target 23:00	90 mg/dL	—	1/50 (2%)	always 100 mg/dL	higher

Kept unchanged in every run

- **Basal rates:** 06:00 0.8 U/h; 10:00 0.7 U/h; 16:00 0.5 U/h; 18:00 0.6 U/h
- **Carb ratios:** 00:00 13 g/U; 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Targets:** 08:00 100 mg/dL
- **AAPS smb_uam:** smb_enabled = true
- **AAPS smb_activation:** smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_high_temp_target_enabled = false
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false

Primary recommendation (first item the model listed)

By area: basal x18, isf x16, aaps_safety_limits x2, dia x2, meal_strategy x1, data_foundation x1, aaps_safety x1, aaps x1, automation_sensitivity_limits x1, uam_aggressiveness x1, UAM SMB delivery x1, automation_sensitivity x1, meal x1, target x1, aaps_uam x1.

- Reduce overnight basal rate x6
- Increase Insulin Sensitivity Factor (ISF) x2
- Reduce overnight basal rate 00:00-05:00 x2
- Begin carbohydrate entries for all meals x1
- Reduce overnight basal (00:00–05:00) x1
- Record all carbohydrate intake consistently x1
- Correct Insulin Sensitivity Factor (ISF) x1
- Reduce overnight basal rate (00:00–05:00) x1
- ... and 34 other titles