

TuneMyBG LLM reproducibility — plain-English summary

Package generated 2026-08-21 covering 2026-08-07 to 2026-08-21 (14 days). Glucose: mean 121.6 mg/dL, time in range 86.0%, below range 5.1%, CV 32.6%. Carb entries logged: 0.

Each run is a brand-new conversation: package attached with Prompt 1, then Prompts 2, 3 and 4 in turn, then the app's correction prompt once if the JSON failed the app's checks. AAPS defaults are from nightscout/AndroidAPS master @ 598e2eb (2026-08-02).

meta-llama/llama-4-maverick — 49 runs (real_llama4maverick_50)

At a glance

- App outcome over 50 completed runs: accepted first time 19; accepted after the repair prompt 25; invalid JSON with no prompt offered 2; schema stop with no prompt offered 0; still rejected after repair 4. Dead ends (no way forward in the app): 2.
- Strict compliance with Prompt 4's own rules (which the app does not fully check): 19/49 (39%) of runs first time.
- Settings changed per run: 0–4 of 23 (most runs: 1).
- 30 settings were left unchanged in every single run; 6 were changed in at least one run, of which 0 were changed in more than half the runs.
- Headline profile recommendation focus: basal in 43/49 (88%), isf in 4/49 (8%), cr in 1/49 (2%), aaps in 1/49 (2%).
- Profile recommendation decision: change ×29, keep ×20.

Changed in at least one run

Frequency is how many of the runs suggested a different value for that setting. Range is the spread of the suggested values across those runs.

Setting	Current	AAPS default	Suggested change in	Suggested values	Direction
Profile / basal 00:00	0.6 U/h	—	19/49 (39%)	0.55–0.65 U/h; most often 0.65 U/h (15 of 19); 2 distinct values	mixed
AAPS / max basal and max IOB limits / max_iob_units	25 U	3.0 U (Simple Mode recalculates this from the profile)	18/49 (37%)	15–20 U; most often 15 U (12 of 18); 2 distinct values	lower
Profile / ISF 00:00	56 mg/dL/U	—	4/49 (8%)	50–55 mg/dL/U; most often 50 mg/dL/U (3 of 4); 2 distinct values	lower
Profile / basal 16:00	0.5 U/h	—	2/49 (4%)	0.55–0.6 U/h; most often 0.55 U/h (1 of 2); 2 distinct values	higher
Profile / basal 06:00	0.8 U/h	—	1/49 (2%)	always 0.75 U/h	lower
Profile / CR 00:00	13 g/U	—	1/49 (2%)	always 12.5 g/U	lower

Kept unchanged in every run

- **Basal rates:** 05:00 0.7 U/h; 10:00 0.7 U/h; 18:00 0.6 U/h; 19:00 0.7 U/h
- **Carb ratios:** 04:00 11.1 g/U; 08:00 13.2 g/U; 16:00 10.5 g/U; 20:00 12.6 g/U
- **Targets:** 00:00 90 mg/dL; 08:00 100 mg/dL; 23:00 90 mg/dL
- **DIA:** 10 h
- **AAPS smb_uam:** smb_enabled = true; uam_enabled = true
- **AAPS smb_activation:** smb_always_enabled = true; smb_with_cob_enabled = false; smb_after_carbs_enabled = false; smb_with_temp_target_enabled = true; smb_with_high_temp_target_enabled = false
- **AAPS smb_delivery:** smb_interval_minutes = 3 min; max_smb_basal_minutes = 15 min; max_uam_smb_basal_minutes = 20 min

- **AAPS safety_limits:** max_basal_u_per_hour = 12 U/h
- **AAPS sensitivity:** autosens_enabled = false; dynamic_sensitivity_enabled = true; autosens_raises_target = false; autosens_lowers_target = false; dynamic_isf_adjustment_factor_percent = 70 %
- **AAPS carb_absorption:** min_5m_carbimpact = 8 mg/dL/5 min; meal_max_absorption_hours = 7 h

Primary recommendation (first item the model listed)

By area: basal ×26, AAPS safety limits ×7, profile ×5, AAPS ×4, meal strategy ×3, Meal Logging ×1, meal logging ×1, Meal Strategy ×1.

- Adjust basal rates ×6
- Verify basal rates ×4
- Adjust max IOB ×4
- Adjust max IOB limit ×3
- Review basal rates ×3
- Verify and adjust basal rates ×2
- Start logging meal data ×2
- Reduce max IOB ×2
- ... and 20 other titles